
## Lite Gait LG150K





Max. Patient Weight	150 kg			
Max. Patient Height Over Ground				
Over Ground	217			
Max. Patient Height w/Treadmill				
Over Treadmill	202			
Yoke Adjustment				
Powered	70			
Manual	N/A			
Overall Unit Height	242			

Base Options	Standard Base <b>76</b>	76/86		
Base Dimensions		76	86	
	Treadmill Width	74	84	
	Door Width	89	99	
Max. Base Footprint 97 x 115				
	Clearance		48	
	Leg Length		115	
Heigh	t (Floor to Top of Leg)	20		
	Caster Diameter	15.2		
Shipping Dimension	ons	107 x 71 x 140 cm		
Shipping Weight		120 kg 125 kg		
Power Requirements: 110V-240V 50/60Hz				

**LiteGait® Includes** Integrated Freedome **INCLUDED Included Harness** 1 Harness BiSym **OPTIONAL** Standard Adult Gaiter Stool **OPTIONAL** 3D Handlebar **INCLUDED** Parellel Grab Bar **INCLUDED OPTIONAL** Q-Straps

Customization Available - Contact Mobility Research Sales Dept. for additional information
Created: 12/12/2018



All dimensions subject to change

Dimensions listed in cm/kg

Price

Treadmill step up height: 16 cm