
Lite Gait LG900





Max. Patient Weight	225 kg			
Max. Patient Height Over Ground				
Over Ground	230			
Max. Patient Height w/Treadmill				
Over Treadmill	215			
Yoke Adjustment				
Powered	90			
Manual	N/A			
Overall Unit Height	255			

Base Options	Standard Base 81	76/81/86		
Base Dimensions		76	81	86
	Treadmill Width	74	79	84
	Door Width	90	95	100
	Max. Base Footprint	97 x 122		
	Clearance	52		
	Leg Length	122		
Height	(Floor to Top of Leg)	20		
	Caster Diameter	15.2		
Shipping Dimensio	ns	107 x 71 x 140 cm		
Shipping Weight		125 kg 130 kg 135 kg		

Power Requirements: 110V-240V 50/60Hz

LiteGait® Includes

Integrated Freedome

INCLUDED

Included Harnesses

Ext. Panel w/43 cm Groin Piece

2 Harness

BiSym

WIRELESS

Small Adult

Gaiter Stool

OPTIONAL

Standard Adult

3D Handlebar

Parellel Grab Bar

INCLUDED

NOT AVAILABLE

Q-Straps

INCLUDED

Customization Available - Contact Mobility Research Sales Dept. for additional information Created: 12/12/2018



All dimensions subject to change Dimensions listed in cm/kg Treadmill step up height: 16 cm

Price