

GaitKeeper GKS22

Model shown with
optional GaitSens 2.0



Walking Surface Width	56
Walking Surface Length	142
Treadmill Width	75
Treadmill Length	185
Footprint	185 L x 75 W
Maximum User Weight	185
Incline Range	0 to 15%
Drive Motor (AC)	0.2 – 16
Reverse	0.2 – 5
Step-Up Height	18
Handrail Length	94
Handrail Height (MAX)	102
AdjustaBars Range	61- 102

Accessories

GaitSens 2.0

Optional

Power Requirements

220V – 230VAC

Regulatory Approvals

50Hz, 10 amp
TUVc, TUVus, UL, CAN/CSA
EMC: CE Mark: Product Specific Standard
IEC60601-1-2:2014-Medical Electrical
Equipment, Electro-Mag Compatibility
LVD: CE Mark: Product Specific Standard
IEC60601-1: 2006+A1: 2012+A12: 2014

Price

All dimensions subject to change.
Dimensions listed in cm/kg/kph.