


# LiteGait® LG700K



Max user weight	150
Max user height overground	217
Max user height over 6" treadmill	202
Powered Stroke adjustment	70
Overall Unit Height	242
Standard base width	76
Other Base Options	76    86
Max treadmill width	74    84
Door width	89    99
Max base footprint	97 x 115
Clearance	48
Leg length	115
Height (floor to top of leg)	20
Caster Diameter	15.2

Accessories		
Integrated FreeDome	Included	
BiSym	Optional	
Gaiter Stool	Optional	
3D Handle Bar	Included	
Parallel Grab Bar	Optional	
Standard Adult Harness	Included	
Q-Straps	Included	

Shipping Dimensions	107 x 71 x 140
Shipping Weight	120   125
Power Requirements	110V-240V 50/60Hz
<b>Price</b>	
[Empty Price Box]	

All dimensions subject to change. Dimensions listed in cm/kg.  
Treadmill step up height: 16 cm

Customizations Available

