


# LiteGait® LG800



Max user weight	200		
Max user height overground	217		
Max user height over 6" treadmill	202		
Powered Stroke adjustment	80		
Overall Unit Height	242		
Standard base width	86		
Other Base Options	76	81	86
Max treadmill width	74	79	84
Door width	89	94	99
Max base footprint	97 x 123		
Clearance	54		
Leg length	123		
Height (floor to top of leg)	22		
Caster Diameter	12.7		

Accessories		
Integrated FreeDome	Included	
BiSym	Not Available	
Gaiter Stool	Optional	
3D Handle Bar	Included	
Parallel Grab Bar	Optional	
Standard Adult Harness	Included	
Q-Straps	Optional	

Shipping Dimensions	107 x 71 x 140
Shipping Weight	125   130   135
Power Requirements	110V-240V 50/60Hz

Price

All dimensions subject to change. Dimensions listed in cm/kg.  
Treadmill step up height: 16 cm

Customizations Available

