


PRODUCT DATA SHEET

LiteGait® LG900



Max user weight	225
Max user height overground	230
Max user height over 6" treadmill	215
Powered Stroke adjustment	90
Overall Unit Height	255
Standard base width (inside)	86
Standard base width (outside)	97
Base options (inside)	76 81 86
Base options (outside)	87 92 97
Max treadmill width	74 79 84
Door width	90 95 100
Max base footprint	97 x 122
Clearance	52
Leg length	122
Height (floor to top of leg)	20
Caster Diameter	15.2

All dimensions subject to change. Dimensions listed in cm/kg.
Treadmill step up height: 16 cm

Customizations Available

Accessories

Integrated FreeDome	Included	
BiSym	Wireless	
Nex Gen Capable	Optional	
Gaiter Stool	Optional	
3D Handle Bar	Included	
Parallel Grab Bar	Not Available	
Standard Adult Harness Ext. panel w/43cm groin piece	Included	
Small Adult Harness	Included	
Q-Straps	Included	

Shipping Dimensions	107 x 71 x 140
Shipping Weight	125 130 135
Power Requirements	110V-240V 50/60Hz

Price



Jordkærvej 29 | 8600 Silkeborg | +45 8680 1807 | info@rehab-partner.dk